




Configure the DNABarcodeCompatibility interface to enable the communication with the R environment

1



When you **first** start the DNABarcodeCompatibility user graphical interface, a window pops up, please select R Settings to set the library paths.



2



On Windows and MacOSx, clicking the R application opens an R console in which to type in **.libPaths()** to list all possible library paths.

Alternatively, on Linux and MacOSx, run R within a terminal and type in **.libPaths()**



3



Copy and paste one-by-one all paths **without quotes** into each empty field of the pop up window

In this pop up window, you can add an empty field by clicking on the "+" icon.

4



Click the save button and then close the pop up explicitly